


Phospho-Pkc Delta-Y311 Rabbit Polyclonal Antibody For Wb, Elisa - Q05655

Item Number: CM0001898



Introduction

Phospho-PKC delta-Y311 Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human, mouse PKC delta. Ideal for studying PKC delta phosphorylation at Tyr311 in apoptosis, cell cycle, and cancer signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC delta-Y311 Rabbit pAb
Remarks/Alias	MAY1; PKCD; ALP53; CVID9; nPKC-delta; Phospho-PKC delta-Y311
Species	Human
GeneID (Human)	5580
GeneID	5580
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around Y311 of rat PKC delta (NP_579841.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q05655
Protein Weight	78kDa
Shipping	Ice bag