

Phospho-Pkc Mu (S916) Rabbit Polyclonal Antibody For Wb And Elisa - Q15139

Item Number: CM0001448

Introduction

Phospho-specific rabbit polyclonal antibody recognizing PRKD1 (PKC mu) phosphorylated at Ser916; suitable for Western blot and ELISA applications. Reacts with human and mouse PRKD1, enabling signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC mu-S916 Rabbit pAb
Remarks/Alias	PKD; PKCM; CHDED; PRKCM; PKC-MU; Phospho-PKC mu-S916
Species	Human
GeneID (Human)	5587
GeneID	5587
Immunogen	A synthetic phosphorylated peptide around S912 of human PKC mu (NP_002733.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15139
Protein Weight	102kDa
Shipping	Ice bag