


Phospho-Pkc (Pan) (Bii Ser660) Rabbit Polyclonal Antibody For Wb, Elisa - P17252 / P05771 / Q05655 / Q02156 / P24723 / Q04759

Item Number: CM0001814



Introduction

High-quality rabbit polyclonal antibody against phosphorylated PKC (pan) (β II Ser660). Validated for WB and ELISA in human, mouse, and rat. Ideal for signal transduction, oncology, and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC (pan) (β II Ser660) Rabbit pAb
Remarks/Alias	PRKCA/PRKCB/PRKCD/PRKCE/PRKCH/PRKCQ; Phospho-PKC (pan) (β II Ser660)
Species	Human
GenelD (Human)	5578/5579/5580/5581/5583/5588
GenelD	5578/5579/5580/5581/5583/5588
Immunogen	Synthetic peptide; A synthetic phosphorylated peptide around S660 of human PKC (pan) β IIPRKCA (NP_002728.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17252/P05771/Q05655/Q02156/P24723/Q04759
Protein Weight	59kDa/67kDa/74- 83kDa
Shipping	Ice bag