

Phospho-PKcalpha/Beta II-T638/641 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P17252 / P05771

Item Number: CM0000518



Introduction

Phospho-PKCalpha/beta II-T638/641 Rabbit pAb detects phosphorylated PKC alpha and beta II in human, mouse, and rat. Suitable for WB, IHC-P, IF/ICC, and ELISA, this antibody supports research in cell signaling, cancer, and cardiovascular biology.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKCalpha/beta II-T638/641 Rabbit pAb
Remarks/Alias	PKC-beta; PKCB; PRKCB1; PRKCB2; Phospho-PKCalpha/beta II-T638/641
Species	Human
GeneID (Human)	5578/5579
GeneID	5578/5579
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T638 of human PKCalpha/beta II PRKCA (NP_002728.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17252/P05771
Protein Weight	76kDa/77kDa
Shipping	Ice bag