


Phospho-Pkcδ/θ-S643/676 Rabbit Pab - Q05655/Q04759

Item Number: CM0001496



Introduction

Phospho-specific rabbit polyclonal antibody targeting PKCδ Ser643 and PKCθ Ser676, validated for WB and ELISA in human and mouse samples. Ideal for studying PKCδ/θ signaling in apoptosis and cancer.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKCδ/θ-S643/676 Rabbit pAb
Remarks/Alias	PRKCD/PRKCQ; Phospho-PKCδ/θ-S643/676
Species	Human
GeneID (Human)	5580/5588
GeneID	5580/5588
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S643 of human PKCδ/θPRKCD (NP_001341605.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q05655/Q04759
Protein Weight	77kDa/80kDa/67kDa/74kDa/81kDa
Shipping	Ice bag