

Phospho-Pkr/Eif2Ak2-T446 Rabbit Mab - P19525

Item Number: CM0001643



Introduction

Rabbit monoclonal antibody targeting phosphorylated T446 of human PKR/EIF2AK2, suitable for WB and ELISA. Ideal for innate immunity, viral infection, and kinase signaling research.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-PKR/EIF2AK2-T446 Rabbit mAb
Remarks/Alias	PKR; PRKR; DYT33; LEUDEN; EIF2AK1; PPP1R83; Phospho-PKR/EIF2AK2-T446
Species	Human
GeneID (Human)	5610
GeneID	5610
Immunogen	Synthetic phosphorylated peptide around T446 of human PKR/EIF2AK2 (P19525)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P19525
Protein Weight	62kDa
Shipping	Ice bag