


Phospho-Plc Gamma 1 (Plcg1)-S1248 Rabbit Monoclonal Antibody For Wb, Elisa - P19174

Item Number: CM0001576



Introduction

Rabbit monoclonal antibody against phosphorylated PLCG1 (S1248), with human and mouse cross-reactivity. Validated for western blot and ELISA. Suitable for studying receptor tyrosine kinase signaling, cell migration, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PLC gamma 1 (PLCG1)-S1248 Rabbit mAb
Remarks/Alias	PLC1; NCKAP3; PLC-II; PLC148; PLCgamma1
Species	Human
GeneID (Human)	5335
GeneID	5335
Immunogen	A synthetic phosphorylated peptide around S1248 of human PLC gamma 1 (PLCG1) (NP_877963.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P19174
Protein Weight	149kDa
Shipping	Ice bag