


Phospho-Plc Gamma 1 (Plcg1)-S1248 Rabbit Pab - P19174

Item Number: CM0000753



Introduction

Rabbit polyclonal antibody detecting phospho-PLCG1 (S1248), validated for WB, IHC-P, ELISA in human, mouse, and rat. Ideal for studying growth factor signaling, actin reorganization, and cell migration.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PLC gamma 1 (PLCG1)-S1248 Rabbit pAb
Remarks/Alias	PLC1; NCKAP3; PLC-II; PLC148; PLCgamma1; Phospho-PLC gamma 1 (PLCG1)-S1248
Species	Human
GeneID (Human)	5335
GeneID	5335
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19174
Protein Weight	149kDa
Shipping	Ice bag