


Phospho-Plc Gamma 1 (Plcg1)-Y783 Rabbit Pab For Wb, Elisa - P19174

Item Number: CM0001231



Introduction

Phospho-PLC gamma 1 (PLCG1)-Y783 Rabbit pAb is a phospho-specific antibody for Western blot and ELISA, targeting phosphorylated Y783 of human PLCG1. Cross-reacts with mouse. Ideal for studying receptor tyrosine kinase signaling, actin reorganization, and cell migration.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PLC gamma 1 (PLCG1)-Y783 Rabbit pAb
Remarks/Alias	PLC1; NCKAP3; PLC-II; PLC148; PLCgamma1; Phospho-PLC gamma 1 (PLCG1)-Y783
Species	Human
GeneID (human)	5335
GeneID	5335
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around Y783 of human PLC gamma 1 (PLCG1) (NP_877963.1?).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P19174
Protein Weight	149kDa
Shipping	Ice bag