


Phospho-Plc Gamma 2 (Plcg2)-Y1217 Rabbit Pab - P16885

Item Number: CM0001783



Introduction

Rabbit polyclonal antibody specific to phosphorylated Y1217 of human PLCG2, validated for WB, IHC-P, ELISA; cross-reacts with human and mouse. Ideal for signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PLC gamma 2 (PLCG2)-Y1217 Rabbit pAb
Remarks/Alias	FCAS3; APLAID; PLC-IV; PLC-gamma-2; Phospho-PLC gamma 2 (PLCG2)-Y1217
Species	Human
GeneID (Human)	5336
GeneID	5336
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16885
Protein Weight	148kDa
Shipping	Ice bag