


Phospho-PLN-S16 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P26678

Item Number: CM0000864



Introduction

Phospho-PLN-S16 Rabbit Monoclonal Antibody (WB, IF-P, IHC-P, ELISA) detects endogenous phospholamban phosphorylated at Ser16 in human, mouse, and rat. Ideal for cardiac calcium homeostasis and cardiomyopathy research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PLN-S16 Rabbit mAb
Remarks/Alias	PLB; CMD1P; CMH18
Species	Human
GeneID (Human)	5350
Immunogen	A synthetic phosphorylated peptide around S16 of human PLN (NP_002658.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P26678
Protein Weight	6 kDa
Shipping	Ice bag