


Phospho-PLN-S16 Rabbit Pab For Wb And Elisa - P26678

Item Number: CM0001651



Introduction

Phospho-PLN-S16 Rabbit pAb specifically detects phosphorylated phospholamban at serine 16 in human, mouse, and rat samples. Ideal for WB and ELISA in cardiac calcium regulation and cardiomyopathy research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PLN-S16 Rabbit pAb
Remarks/Alias	PLB; CMD1P; CMH18; Phospho-PLN-S16
Species	Human
GeneID (Human)	5350
GeneID	5350
Immunogen	A synthetic phosphorylated peptide around S16 of human PLN (NP_002658.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26678
Protein Weight	6kDa
Shipping	Ice bag