


## Phospho-PLN-S16/T17 Rabbit Polyclonal Antibody For Wb, Elisa - P26678

Item Number: CM0001541



### Introduction

Phospho-PLN-S16/T17 rabbit polyclonal antibody for Western blot and ELISA, reacts with human, mouse, rat. Detects phosphorylated phospholamban, a key modulator of SERCA2 activity in cardiac muscle sarcoplasmic reticulum.

[Learn More](#)

| Parameter        | Value                                                                           |
|------------------|---------------------------------------------------------------------------------|
| Product Name     | Phospho-PLN-S16/T17 Rabbit pAb                                                  |
| Remarks/Alias    | PLB; CMD1P; CMH18; Phospho-PLN-S16/T17                                          |
| Species          | Human                                                                           |
| Gene ID (Human)  | 5350                                                                            |
| Gene ID          | 5350                                                                            |
| Immunogen        | A synthetic phosphorylated peptide around S16 & T17 of human PLN (NP_002658.1). |
| Source           | Rabbit                                                                          |
| Category         | Phosphorylated Antibodies                                                       |
| Application      | WB, ELISA                                                                       |
| Cross Reactivity | Human, Mouse, Rat                                                               |
| SWISS            | P26678                                                                          |
| Protein Weight   | 6kDa                                                                            |
| Shipping         | Ice bag                                                                         |