


Phospho-PLN-T17 Rabbit Pab - P26678

Item Number: CM0001871



Introduction

High-specificity Phospho-PLN-T17 rabbit polyclonal antibody for WB, IHC-P, ELISA. Detects human, mouse, rat phospholamban phosphorylated at Thr17. Ideal for cardiac calcium signaling and cardiomyopathy research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PLN-T17 Rabbit pAb
Remarks/Alias	PLB; CMD1P; CMH18; Phospho-PLN-T17
Species	Human
GeneID (human)	5350
GeneID	5350
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T17 of human PLN (NP_002658.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26678
Protein Weight	6kDa
Shipping	Ice bag