


Phospho-Polr2A Ctd-S2 Rabbit Mab - P24928

Item Number: CM0001971



Introduction

Rabbit monoclonal antibody specific for phosphorylated Ser2 of POLR2A CTD. Validated for WB, IF/ICC, IHC-P, ChIP, ChIP-seq, CUT&Tag, ELISA; reacts with human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	Phospho-POLR2A CTD-S2 Rabbit mAb
Remarks/Alias	RPB1; RPO2; POLR2; POLRA; RPBh1; RPOL2; NEDHIB; RpIIL5; hsRPB1; hRPB220; Phospho-POLR2A CTD-S2
Species	Human
GeneID (Human)	5430
GeneID	5430
Immunogen	Synthetic Peptide (A phospho specific peptide corresponding to residues surrounding S2 of human POLR2A CTD repeat YSPTSPS)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, IHC-P, ChIP, ChIP-seq, CUT&Tag, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24928
Protein Weight	217 kDa
Shipping	Ice bag