


Phospho-Polr2A Ctd- S2/S5 Rabbit Pab - P24928

Item Number: CM0001491



Introduction

Phospho-POLR2A CTD (S2/S5) rabbit polyclonal antibody validated for WB and ELISA. Targets human POLR2A CTD repeat YSPTSPS phosphorylation at Ser2/Ser5. Ideal for transcription research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-POLR2A CTD- S2/S5 Rabbit pAb
Remarks/Alias	RPB1; RPO2; POLR2; POLRA; RPBh1; RPOL2; NEDHIB; RplIL5; hsRPB1; hRPB220; Phospho-POLR2A CTD- S2/S5
Species	Human
GeneID (Human)	5430
GeneID	5430
Immunogen	Synthetic peptide A phospho specific peptide corresponding to residues surrounding S2/S5 of human POLR2A CTD repeat YSPTSPS.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P24928
Protein Weight	217kDa
Shipping	Ice bag