


Phospho-Polr2A Ctd-S5 Rabbit Mab For Wb, Ihc-P, If/Icc, Elisa - P24928

Item Number: CM0000436



Introduction

Phospho-specific rabbit monoclonal antibody against POLR2A CTD Ser5, validated for WB, IHC-P, IF/ICC, and ELISA in human, mouse, and rat. Ideal for studying RNA polymerase II transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-POLR2A CTD-S5 Rabbit mAb
Remarks/Alias	RPB1; RPO2; POLR2; POLRA; RPBh1; RPOL2; NEDHIB; RpIIL5; hsRPB1; hRPB220; Phospho-POLR2A CTD-S5
Species	Human
Gene ID (Human)	5430
Gene ID	5430
Immunogen	Synthetic Peptide A phospho specific peptide corresponding to residues surrounding S5 of human POLR2A CTD (P24928).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24928
Protein Weight	217kDa
Shipping	Ice bag