


Phospho-Polr2A Ctd-T4 Rabbit Polyclonal Antibody For Wb, Ip, Elisa, Chip - P24928

Item Number: CM0000458



Introduction

High-quality Phospho-POLR2A CTD-T4 Rabbit polyclonal antibody validated for WB, IP, ELISA, ChIP; detects human, mouse, rat POLR2A phosphorylated at T4 of CTD heptapeptide repeats. Ideal for transcription studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-POLR2A CTD-T4 Rabbit pAb
Remarks/Alias	RPB1; RPO2; POLR2; POLRA; RPBh1; RPOL2; NEDHIB; RplILS; hsRPB1; hRPB220; Phospho-POLR2A CTD-T4
Species	Human
Gene ID (Human)	5430
Gene ID	5430
Immunogen	Synthetic Peptide A phospho specific peptide corresponding to residues surrounding T4 of human POLR2A CTD repeat YSPTSPS.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	P24928
Protein Weight	217kDa
Shipping	Ice bag