


Phospho-Pp2A Catalytic A (Y307) Rabbit Monoclonal Antibody For Wb, Elisa - P67775

Item Number: CM0001595



Introduction

Rabbit monoclonal antibody specific to phospho-PP2A catalytic α (Tyr307) for Western blot and ELISA. Recognizes human and mouse PPP2CA; predicted molecular weight 36 kDa. SWISS: P67775.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PP2A Catalytic α -Y307 Rabbit mAb
Remarks/Alias	RP-C; PP2Ac; PP2CA; NEDLBA; PP2Calpha; Phospho-PP2A Catalytic α -Y307
Species	Human
Gene ID (Human)	5515
Immunogen	Synthetic Peptide - A synthetic phosphorylated peptide around Y307 of human PP2A Catalytic α (P67775)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P67775
Protein Weight	36kDa
Shipping	Ice bag