


Phospho-Ppig-S376 Rabbit Pab For Wb, Elisa - Q13427

Item Number: CM0000953



Introduction

Phospho-specific rabbit polyclonal antibody targeting PPIG phosphorylated at S376, validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for studying splicing regulation and phospho-signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PPIG-S376 Rabbit pAb
Remarks/Alias	CYP; SRCyp; SCAF10; CARS-Cyp; Phospho-PPIG-S376
Species	Human
Gene ID	9360
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S376 of human PPIG (NP_004783.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13427
Protein Weight	89kDa
Shipping	Ice bag