


Phospho-Ppp1R12A-S668 Rabbit Polyclonal Antibody For Wb, Elisa - O14974

Item Number: CM0000885



Introduction

Phospho-PPP1R12A (MYPT1) S668 rabbit polyclonal antibody validated for WB and ELISA. Detects PPP1R12A phosphorylated at serine 668 in human, mouse, and rat. Ideal for myosin phosphatase signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PPP1R12A-S668 Rabbit pAb
Remarks/Alias	MBS; GUBS; M130; MYPT1; Phospho-PPP1R12A-S668
Species	Human
Gene ID	4659
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14974
Protein Weight	115kDa
Shipping	Ice bag