


Phospho-Ppp1R12A/Ppp1R12B/Ppp1R12C (Thr696) Rabbit Polyclonal Antibody For Wb, Elisa - O14974 / O60237 / Q9Bzl4

Item Number: CM0001338



Introduction

Rabbit polyclonal antibody specific to phosphorylated Thr696 of human PPP1R12A/PPP1R12B/PPP1R12C myosin phosphatase regulatory subunits. Validated for Western blot and ELISA in human samples. Ideal for cell signaling and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PPP1R12A/PPP1R12B/PPP1R12C-T696 Rabbit pAb
Species	Human
GeneID (Human)	4659/4660/54776
GeneID	4659/4660/54776
Immunogen	Synthetic phosphorylated peptide around T696 of human PPP1R12A/PPP1R12B/PPP1R12C (based on PPP1R12A NP_002471.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14974 / O60237 / Q9BZL4
Shipping	Ice bag