


Phospho-Prkcq-S676 Rabbit Pab - Q04759

Item Number: CM0001811



Introduction

Phospho-specific rabbit polyclonal antibody targeting PRKCQ (S676). Validated for WB and ELISA detection of human phosphorylated PRKCQ, a key kinase in T-cell signaling, inflammation, and platelet function.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PRKCQ-S676 Rabbit pAb
Remarks/Alias	PRKCT; nPKC-theta; Phospho-PRKCQ-S676
Species	Human
Gene ID (Human)	5588
Gene ID	5588
Immunogen	Synthetic peptide corresponding to residues surrounding S676 of human PRKCQ (phospho-specific)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q04759
Protein Weight	82kDa
Shipping	Ice bag