


Phospho-Prkcq-T538 Rabbit Polyclonal Antibody For Wb, Elisa - Q04759

Item Number: CM0001688



Introduction

Phospho-specific rabbit polyclonal antibody detecting human PRKCQ phosphorylated at Thr538. Validated for WB and ELISA, this antibody targets protein kinase C theta, a key kinase in T-cell receptor signaling, NF-kB activation, and immune response.

[Learn More](#)

Core Product Specifications	
Product Name	Phospho-PRKCQ-T538 Rabbit pAb
Remarks/Alias	PRKCT; nPKC-theta; Phospho-PRKCQ-T538
Species	Human
Gene ID (Human)	5588
Gene ID	5588
Immunogen	A synthetic phosphorylated peptide around T538 of human PRKCQ (NP_006248.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q04759
Protein Weight	82kDa
Shipping	Ice bag