


Phospho-Prmt4/Carm1 (Ser228) Rabbit Polyclonal Antibody For Wb And Elisa - Q86X55

Item Number: CM0001277



Introduction

Rabbit polyclonal antibody recognizing phospho-CARM1 (Ser228) for WB and ELISA. Detects human CARM1 (~66 kDa). Useful for studying transcriptional coactivation, RNA processing, and cancer signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PRMT4/CARM1-S228 Rabbit pAb
Remarks/Alias	PRMT4; Phospho-PRMT4/CARM1-S228
Species	Human
GeneID (Human)	10498
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q86X55
Protein Weight	66kDa
Shipping	Ice bag