

Phospho-Psd93/Chapsyn-110/Dlg2-Y340 Rabbit Pab For Wb And Elisa - Q15700

Item Number: CM0000901



Introduction

Phospho-PSD93/DLG2-Y340 rabbit polyclonal antibody for WB, ELISA. Validated in human, mouse, rat. Detects postsynaptic scaffold protein chapsyn-110/PSD-93 involved in NMDA receptor signaling and chronic pain perception.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PSD93/chapsyn-110/DLG2-Y340 Rabbit pAb
Remarks/Alias	PSD93; PSD-93; PPP1R58; chapsyn-110; Phospho-PSD93/chapsyn-110/DLG2-Y340
Species	Human
Gene ID (Human)	1740
Gene ID	1740
Immunogen	A synthetic phosphorylated peptide around Y340 of human PSD93/chapsyn-110/DLG2 (NP_001355.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15700
Protein Weight	98kDa
Shipping	Ice bag