


Phospho-Rac1/2/3-S71 Rabbit Pab For Wb, Elisa - P63000/P15153/P60763

Item Number: CM0001411



Introduction

Phospho-Rac1/2/3-S71 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects Rac1/2/3 phosphorylated at S71 in human and rat, 21 kDa. Essential for studying cell migration, adhesion, and NADPH oxidase activity.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Rac1/2/3-S71 Rabbit pAb
Remarks/Alias	MIG5; MRD48; Rac-1; TC-25; p21-Rac1; Phospho-Rac1/2/3-S71
Species	Human
GeneID (Human)	5879/5880/5881
GeneID	5879/5880/5881
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S71 of human RAC1 (NP_008839.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P63000/P15153/P60763
Protein Weight	21kDa
Shipping	Ice bag