

Phospho-Raf1-S338 Rabbit Pab For Wb, Ihc-P, Elisa - P04049

Item Number: CM0000994



Introduction

Rabbit polyclonal antibody specific for phosphorylated Raf1 at Ser338 (Swiss-Prot: P04049). Validated for western blot, IHC-P, and ELISA. Reacts with human, mouse, and rat Raf1. Suitable for MAPK signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Raf1-S338 Rabbit pAb
Remarks/Alias	NS5; CRAF; Raf-1; c-Raf; CMD1NN; Phospho-Raf1-S338
Species	Human
GeneID (Human)	5894
GeneID	5894
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S338 of human Phospho-Raf1-S338 (NP_002871.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04049
Protein Weight	73kDa
Shipping	Ice bag