

Phospho-Raptor-Ser792 Rabbit Monoclonal Antibody For Wb, Elisa - Q8N122

Item Number: CM0000746



Introduction

Phospho-Raptor-Ser792 Rabbit Monoclonal Antibody validated for WB and ELISA. Detects human RPTOR, a regulatory component of mTORC1 complex essential for nutrient sensing and cell growth signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Raptor-Ser792 Rabbit mAb
Remarks/Alias	KOG1; Mip1
Species	Human
Gene ID (Human)	57521
Gene ID	57521
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8N122
Protein Weight	149 kDa
Shipping	Ice bag