


Phospho-Rb-S811 Rabbit Pab For Wb, Elisa - P06400

Item Number: CM0000473



Introduction

Phospho-specific rabbit polyclonal antibody targeting human RB phosphorylated at serine 811, validated for Western blot and ELISA. Ideal for research on cell cycle regulation, tumor suppression, and E2F signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-RB-S811 Rabbit pAb
Remarks/Alias	RB; pRb; OSRC; pp110; p105-Rb; PPP1R130; p110-RB1; Phospho-RB-S811
Species	Human
GeneID (Human)	5925
GeneID	5925
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S811 of human RB (NP_000312.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P06400
Protein Weight	106kDa
Shipping	Ice bag