


Phospho-Relb-S573 Rabbit Polyclonal Antibody For Wb And Elisa - Q01201

Item Number: CM0000618



Introduction

Rabbit polyclonal antibody targeting phosphorylated RELB at S573. Validated for WB and ELISA; cross-reacts with mouse. Useful for investigating NF-kB signaling, immune cell function, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	Phospho-RELB-S573 Rabbit pAb
Remarks/Alias	IREL; I-REL; IMD53; REL-B; Phospho-RELB-S573
Species	Human
GenelD	5971 (Human)
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q01201
Protein Weight	62kDa
Shipping	Ice bag