


Phospho-Rpa2-S4/S8 Rabbit Pab - P15927

Item Number: CM0001088



Introduction

Phospho-specific rabbit polyclonal antibody targeting RPA2 phosphorylated at S4/S8. Ideal for WB and ELISA applications. Detects human RPA2 and cross-reacts with rat. Suitable for DNA replication and damage response studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-RPA2-S4/S8 Rabbit pAb
Remarks/Alias	REPA2; RPA32; RP-A p32; RP-A p34; Phospho-RPA2-S4/S8
Species	Human
Gene ID (Human)	6118
Gene ID	6118
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S4 & S8 of human RPA2 (NP_002937.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P15927
Protein Weight	29kDa
Shipping	Ice bag