


Phospho-Rpb1 Ctd-S2/S5 Rabbit Mab - P24928

Item Number: CM0001654



Introduction

Phospho-specific rabbit monoclonal antibody targeting RNA polymerase II CTD Ser2/Ser5 phosphorylation. Validated for WB, IP, ChIP, CUT&Tag, and ELISA in human samples. Ideal for transcription research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Rpb1 CTD-S2/S5 Rabbit mAb
Remarks/Alias	RPB1; RPO2; POLR2; POLRA; RPBh1; RPOL2; NEDHIB; RplIL5; hsRPB1; hRPB220; Phospho-POLR2A-S2/S5
Species	Human
GeneID (human)	5430
GeneID	5430
Immunogen	A phospho specific peptide corresponding to residues surrounding S2/S5 of human POLR2A
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ChIP, CUT&Tag, ELISA
Cross Reactivity	Human
SWISS	P24928
Protein Weight	217kDa
Shipping	Ice bag