


Phospho-Rps6Ka5-T581 Rabbit Pab - O75582

Item Number: CM0001492



Introduction

Phospho-RPS6KA5-T581 rabbit polyclonal antibody specifically detects the phosphorylated T581 site of human RPS6KA5 (MSK1). Validated for WB and ELISA, with cross-reactivity to rat. Ideal for studying stress and mitogen-activated signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-RPS6KA5-T581 Rabbit pAb
Remarks/Alias	MSK1; RLPK; MSPK1; Phospho-RPS6KA5-T581
Species	Human
Gene ID (Human)	9252
Gene ID	9252
Immunogen	Synthetic peptide, A synthetic phosphorylated peptide around T581 of human RPS6KA5 (NP_004746.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	O75582
Protein Weight	90kDa
Shipping	Ice bag