

Phospho-S6 Ribosomal Protein (Rps6)-S235/236 Rabbit Mab - P62753

Item Number: CM0000831



Introduction

Phospho-RPS6 (S235/236) rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Detects 29 kDa phosphorylated ribosomal protein S6 in human, mouse, rat. Ideal for studying cell growth, proliferation, and translational control.

[Learn More](#)

Parameter	Value
Product Name	Phospho-S6 Ribosomal Protein (RPS6)-S235/236 Rabbit mAb
Remarks/Alias	S6; eS6; Phospho-S6 Ribosomal Protein (RPS6)-S235/236
Species	Human
Gene ID (Human)	6194
Gene ID	6194
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S235 & S236 of human S6 Ribosomal Protein (RPS6) (NP_001001.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62753
Protein Weight	29 kDa
Shipping	Ice bag