

Phospho-S6 Ribosomal Protein (Rps6)-S240/244 Rabbit Mab For Wb, Ihc-P, Elisa - P62753

Item Number: CM0000522



Introduction

Phospho-specific rabbit monoclonal antibody against RPS6 (Ser240/244), validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, and rat RPS6. Ideal for mTOR signaling and ribosome biogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-S6 Ribosomal Protein (RPS6)-S240/244 Rabbit mAb
Remarks/Alias	S6; eS6; Phospho-S6 Ribosomal Protein (RPS6)-S240/244
Species	Human
GeneID (human)	6194
GeneID	6194
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S240 & S244 of human Phospho-S6 Ribosomal Protein (RPS6) (NP_001001.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62753
Protein Weight	29kDa
Shipping	Ice bag