

Phospho-(Ser/Thr) Atm/Atr Substrate Rabbit Pab

Item Number: CM0000639

Introduction

Rabbit polyclonal antibody targeting phosphorylated Ser/Thr in ATM/ATR substrates. Validated for Western blot and ELISA in human, mouse, and rat samples. Ideal for studying DNA damage signaling and cell cycle checkpoints.

[Learn More](#)

Parameter	Value
Product Name	Phospho-(Ser/Thr) ATM/ATR Substrate Rabbit pAb
Species	Human
Gene ID (Human)	pan p-S/T
Immunogen	A synthetic peptide corresponding to a sequence containing phosphorylated S & T.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Shipping	Ice bag