

Phospho-Sf3B1-T313 Rabbit Pab For Wb, Elisa - 075533

Item Number: CM0000800



Introduction

Phospho-SF3B1 (T313) rabbit polyclonal antibody validated for WB and ELISA. Cross-reacts with human, mouse, and rat SF3B1 (~146 kDa). Ideal for spliceosome and pre-mRNA splicing research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-SF3B1-T313 Rabbit pAb
Remarks/Alias	MDS; PRP10; Hsh155; PRPF10; SAP155; SF3b155; Phospho-SF3B1-T313
Species	Human
GeneID (Human)	23451
GeneID	23451
Immunogen	Synthetic Peptide – A synthetic phosphorylated peptide around T313 of human SF3B1 (NP_036565.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075533
Protein Weight	146kDa
Shipping	Ice bag