


Phospho-Shc1-Y427 Rabbit Pab - P29353

Item Number: CM0001636



Introduction

Phospho-SHC1-Y427 rabbit polyclonal antibody for WB and ELISA, reacts with human and mouse SHC1. Specific detection of Y427 phosphorylation for studying SHC1 signaling, angiogenesis, and oxidative stress response.

[Learn More](#)

Parameter	Value
Product Name	Phospho-SHC1-Y427 Rabbit pAb
Remarks/Alias	SHC; SHCA; Phospho-SHC1-Y427
Species	Human
GeneID (Human)	6464
GeneID	6464
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y427 of human SHC1 (NP_892113.4)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P29353
Protein Weight	63kDa
Shipping	Ice bag