

Phospho-Ship1-Y1020 Rabbit Pab For Wb, Elisa - Q92835

Item Number: CM0001928



Introduction

Phospho-SHIP1-Y1020 Rabbit Polyclonal Antibody detects phosphorylated human SHIP1 at Tyrosine 1020 for Western blot and ELISA. Ideal for studying negative regulation of PI3K pathways in immune and hematopoietic cells.

[Learn More](#)

Parameter	Value
Product Name	Phospho-SHIP1-Y1020 Rabbit pAb
Remarks/Alias	SHIP; SHIP1; SHIP-1; hp51CN; SIP-145; p150Ship; Phospho-SHIP1-Y1020
Species	Human
Gene ID (Human)	3635
Gene ID	3635
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y1020 of human SHIP1 (NP_001017915.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q92835
Protein Weight	133kDa
Shipping	Ice bag