


Phospho-Shp2-Y542 Rabbit Pab - Q06124

Item Number: CM0001522



Introduction

Phospho-SHP2-Y542 Rabbit pAb detects SHP2 phosphorylated at tyrosine 542 in human, mouse, rat samples. Validated for Western blot (WB) and ELISA. Essential for MAPK signaling and cancer research. - Q06124

[Learn More](#)

Parameter	Specification
Product Name	Phospho-SHP2-Y542 Rabbit pAb
Remarks/Alias	CFC; NS1; JMML; SHP2; BPTP3; PTP2C; METCDS; PTP-1D; SH-PTP2; SH-PTP3; Phospho-SHP2-Y542
Species	Human
Gene ID (Human)	5781
Gene ID	5781
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06124
Protein Weight	68kDa
Shipping	Ice bag