

Phospho-Sirt1-S47 Rabbit Mab For Western Blot And Elisa - Q96Eb6

Item Number: CM0001125



Introduction

Phospho-SIRT1 (Ser47) Rabbit Monoclonal Antibody validated for WB and ELISA. Detects human SIRT1 at ~82 kDa. Ideal for studying NAD-dependent deacetylation, metabolism, apoptosis, and DNA damage response.

[Learn More](#)

Parameter	Value
Product Name	Phospho-SIRT1-S47 Rabbit mAb
Remarks/Alias	SIR2; SIR2L1; SIR2alpha; Phospho-SIRT1-S47
Species	Human
Gene ID (Human)	23411
Gene ID	23411
Immunogen	A synthetic phosphorylated peptide around S47 of human SIRT1 (Q96EB6)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96EB6
Protein Weight	82kDa
Shipping	Ice bag