


Phospho-Smad1-S206 Rabbit Pab For Wb, Elisa - Q15797

Item Number: CM0001085



Introduction

Phospho-Smad1 (Ser206) rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat Smad1 phosphorylated at serine 206. Suitable for BMP signaling pathway research and phosphorylation studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Smad1-S206 Rabbit pAb
Remarks/Alias	BSP1; JV41; BSP-1; JV4-1; MADH1; MADR1; Phospho-Smad1-S206
Species	Human
Gene ID (Human)	4086
Gene ID	4086
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S206 of human Smad1 (NP_005891.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15797
Protein Weight	52kDa
Shipping	Ice bag