

Phospho-Smad2 (S245/S250/S255) Rabbit Mab For Wb And Elisa - Q15796

Item Number: CM0001588



Introduction

Phospho-Smad2-S245/S250/S255 rabbit monoclonal antibody validated for Western blot and ELISA in human samples. Recognizes phosphorylated SMAD2, a key TGF-beta signaling transducer. Ideal for studying signal transduction, tumor suppression, and developmental biology.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Smad2-S245/S250/S255 Rabbit mAb
Remarks/Alias	JV18; LDS6; CHTD8; MADH2; MADR2; JV18-1; hMAD-2; hSMAD2; Phospho-Smad2-S245/S250/S255
Species	Human
GeneID (Human)	4087
GeneID	4087
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S245 & S250 & S255 of human Smad2 Q15796
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15796
Protein Weight	52 kDa
Shipping	Ice bag