

Phospho-Smad2-S465/467 + Smad3-S423/425 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q15796 / P84022

Item Number: CM0001236



Introduction

Rabbit polyclonal antibody against dual phosphorylated Smad2 (S465/467) and Smad3 (S423/425). Validated for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for TGF-beta signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Smad2-S465/467 + Smad3-S423/425 Rabbit pAb
Remarks/Alias	HSPC193; Hst17436; JV15-2; LDS1C; LDS3; MADH3; Phospho-Smad2-S465/467 + Smad3-S423/425
Species	Human
Gene ID (Human)	4087/4088
Gene ID	4087/4088
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S465 & S467 of human Phospho-Smad2-S465/467 + Smad3-S423/425SMAD2 (NP_005892.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15796/P84022
Protein Weight	48kDa/52kDa/25kDa/35kDa/43kDa
Shipping	Ice bag