


Phospho-Smad2-S465/S467 Rabbit Pab For Wb, Elisa - Q15796

Item Number: CM0000494



Introduction

Phospho-Smad2-S465/S467 rabbit polyclonal antibody for WB and ELISA. Targets TGF-beta-activated Smad2 in human and mouse samples (~52 kDa). Ideal for signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Smad2-S465/S467 Rabbit pAb
Remarks/Alias	JV18; LDS6; CHTD8; MADH2; MADR2; JV18-1; hMAD-2; hSMAD2; Phospho-Smad2-S465/S467
Species	Human
Gene ID (Human)	4087
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S465 & S467 of human Smad2 (NP_005892.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15796
Protein Weight	52 kDa
Shipping	Ice bag