


Phospho-Smad2-T220 Rabbit Pab - Q15796

Item Number: CM0000743



Introduction

Phospho-Smad2-T220 Rabbit pAb is a rabbit polyclonal antibody that detects endogenous levels of SMAD2 only when phosphorylated at threonine 220. Validated for WB, IHC-P, ELISA, it reacts with human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Smad2-T220 Rabbit pAb
Remarks/Alias	JV18; LDS6; CHTD8; MADH2; MADR2; JV18-1; hMAD-2; hSMAD2; Phospho-Smad2-T220
Species	Human
GeneID (human)	4087
GeneID	4087
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around T220 of human Smad2 (NP_005892.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15796
Protein Weight	52kDa
Shipping	Ice bag