


Phospho-Smc1-S957 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q14683

Item Number: CM0001677



Introduction

Phospho-specific SMC1 pS957 rabbit monoclonal antibody for WB, IHC-P, ELISA. Reacts with human, mouse, rat SMC1A. Ideal for studying chromosome cohesion and DNA repair pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-SMC1-S957 Rabbit mAb
Remarks/Alias	SMC1; SMCB; CDLS2; DEE85; SB1.8; EIEE85; SMC1L1; DXS423E; SMC1alpha; Phospho-SMC1-S957
Species	Human
GeneID	8243
Immunogen	A synthetic phosphorylated peptide around S957 of human SMC1 (Q14683)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14683
Protein Weight	143kDa
Shipping	Ice bag