

Phospho-Smo-S616 Rabbit Pab - Q99835

Item Number: CM0000448



Introduction

High-quality rabbit polyclonal antibody against phosphorylated Smoothed (SMO) at Ser616. Validated for WB, DB, and ELISA in human samples, ideal for hedgehog signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-SMO-S616 Rabbit pAb
Remarks/Alias	Gx; CRJS; PHLS; SMOH; FZD11; Phospho-SMO-S616
Species	Human
GeneID (Human)	6608
GeneID	6608
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human
SWISS	Q99835
Protein Weight	86kDa
Shipping	Ice bag