

Phospho-Sp1 (T739) Rabbit Polyclonal Antibody For Western Blotting - P08047

Item Number: CM0000921



Introduction

Rabbit polyclonal antibody specific for phosphorylated Sp1 (Thr739). Validated for Western blot in human, mouse, and rat. Sp1 is a key transcription factor regulating cell growth, apoptosis, and circadian rhythm.

[Learn More](#)

Parameter	Value
Product Name	Phospho-SP1-T739 Rabbit pAb
Remarks/Alias	SP1; Phospho-SP1-T739
Species	Human
Gene ID	6667
Immunogen	A phospho specific peptide corresponding to residues surrounding T739 of human SP1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human, Mouse, Rat
SWISS	P08047
Protein Weight	81kDa
Shipping	Ice bag